

Get Breast Milk Home from Anywhere in the World

portal.milkstork.com/benesch




How It Works



PUMP


Pump according to your regular nursing schedule and refrigerate or freeze your breast milk.



PACK

Pack your milk into the cooler provided.

If shipping, your cooler will be pre-labeled with a FedEx Priority Overnight shipping label.



SHIP





TOTE/CHECK


Ship your milk from an authorized FedEx shipping location or schedule a FedEx pick-up with your hotel or destination office.

Tote your cooler or check it as luggage using the Milk Stork travel bags provided.

Let Milk Stork Do the Heavy Lifting




Simple Order Process & Travel Logistics



Pre-Labeled Ready-To-Go Refrigerated & Frozen Kits



Expedited Shipping on All Deliveries



Proactive Shipment
Monitoring

Need more info? Visit <u>support.milkstork.com</u>